



Background

- Tumor expression of PD-L1 (B7H1) is associated with a response to both anti-PD-L1 and anti-PD-1 therapy.
- Agents that inhibit the PD-1/PD-L1 interaction have been approved in melanoma and are in clinical trials for many other indications, including NSCLC.
- OncoPlexDx offers a quantitative PD-L1 assay to identify PD-L1 positive tumors that might respond to anti-PD-1/PD-L1 therapy.




Figure 1. Immunological Synapse between T-cell and Tumor

- Besides PD-1/PD-L1 therapy, other immuno-oncology drugs have been approved. Additionally, other protein pairs (Figure 1) that promote tumor survival through immune-system evasion are being evaluated as druggable targets.
- OncoPlexDx is currently building an immuno-oncology panel to screen FFPE patient biopsies for expression of druggable immuno targets, which can be assayed simultaneously with our targeted therapy panel.

Methods




Figure 2. Liquid Tissue® - SRM workflow for protein analysis from FFPE tissue.

- Optimal quantitative peptides were identified using recombinant PD-L1, and the assay was validated and characterized using 14 formalin fixed cancer cell lines and FFPE clinical tissue.

Results




Figure 3. Correlation between PD-L1 protein and PD-L1 mRNA. Fourteen cancer cell lines were screened for PD-L1 protein expression, and compared to mRNA expression (Barretina et al., Nature 483:603-607). PD-L1 protein expression was detected in 7 cell lines with a moderate correlation with PD-L1 mRNA levels, suggesting that analysis of mRNA is not representative of protein levels.

Results (continued)

Incidence of PD-L1 Expression Increases With Cancer Stage

The PD-L1 assay was run on archival FFPE sections from 9 normal lung tissues, and 31 early staged (stage 1 and stage 2) and 8 advanced staged (Stage 3) NSCLC tumors (Table 1).

Table 1. Sample classification and PD-L1 protein expression

	Normal	Early Stage (Stage 1 and 2)		Advanced Stage (Stage 3)	
		Squamous cell carcinoma	Adeno carcinoma	Squamous cell carcinoma	Adeno carcinoma
Sample #	9	25	6	5	3
PD-L1 +	0	6	0	2	1
%	0 %	24 %	0 %	40 %	33 %




Figure 4. PD-L1 protein expression in normal and NSCLC samples: Advanced NSCLC patients are more likely to be PD-L1 positive compared to early stage NSCLC patients.

Abundance Level of PD-L1 Protein in Archival NSCLC tumors





Figure 4. Distribution of PD-L1 protein expression. PD-L1 was detected in ~23% of NSCLC samples. The mean(±SD) of the PD-L1 positive samples is demarcated with black bars.

Identification of PD-L1 Positive Patient Tumors

Tumor Characteristics	PDL1 amol/µg
High grade transitional cell carcinoma of bladder	318
Poorly differentiated lung adenocarcinoma	343
Metastatic adenocarcinoma of lung	351
Malignant neoplasm of lung	184
Stomach adenocarcinoma	1025
Metastatic moderately differentiated adenocarcinoma of primary gastric origin	349



Immuno-Oncology Panel in Development

Immuno-Oncology (IO) Panel	Tumor Infiltrating Lymphocytes (TIL) Panel
PD-L1, PD-L2, B7.1, B7.2, B7H2, B7H3, B7H4, CD70, CD40, CD137L, OX40L, IDO1, GAL9, HVEM	PD-1, CTLA4, CD137, OX40, CD27, CD40L, VISTA, BTLA, CD160, LAG3, TIM3, CD28, CD3, CD8




Figure 5. Diagnostic panel to profile Tumor-TIL interaction. The panel can be used to make informed decisions about the immuno-oncology based combination therapy. Assays have been developed to proteins highlighted in red. Additional protein targets are currently being evaluated

Targeted Therapy + Immuno-Therapy

- Molecular diagnostics-based combination therapy (targeted therapy & immunotherapy) is used to reduce toxicity and increase effectiveness.
- Only targeted proteomics (LT-SRM) can simultaneously quantify multiple proteins in FFPE patient biopsies with limited tissue.
- Additional available assays that can help pair therapies are shown in table below.

ALK	ROS1	RET	EGFR	HER2	HER3	HER4	MET	IGF1R	RON
RRM1	hENT1	ERCC1	XRCC1	TOPO1	TOPO2A	TS	TUBB3	SPARC	AR
FRalpha	AXL	FGFR1	FGFR2	FGFR3	BRAF	INSR	PSMA	MSLN	IDO1

Discussion

- We have developed a quantitative mass spectrometry based PD-L1 assay for clinical samples. This assay can be multiplexed with other immuno targets to better understand tumor-T-cell biology.
- Targeted proteomics can quantitate multiple proteins from both the tumor and the tumor infiltrating lymphocytes using only a few sections of FFPE tissue.
- The OncoPlexDx platform allows for quantitation of multiple actionable protein targets for optimal therapy selection based on tumor biology.