

Sarit Schwartz,¹ Shankar Sellappan,¹ Steven W Mamus,² Daniel VT Catenacci,³ Adele Blackler,¹ Fabiola Cecchi,¹ Todd Hembrough¹

¹NantOmics., Rockville, MD ; ²Cancer Center of Sarasota-Manatee, Sarasota, FL; ³University of Chicago, Department of Medicine, Section of Hematology and Oncology, Chicago, IL

Background

- Successful, personalized cancer treatment requires the integration of both genomics and proteomics profiling to identify which therapies will likely be most beneficial for an individual patient, based on the biomarker fingerprint for each patient's unique tumor.
- Patients with HER2-positive breast and gastric cancers have responded remarkably well to anti-HER2 therapy, prompting ongoing trials of such therapies in other indications.
- NantOmics multiplexed mass spectrometry (MS)-based proteomic analysis objectively quantifies HER2 and other actionable biomarkers from two FFPE tissue sections. High HER2 protein levels (≥ 740 and ≥ 750 amol/ug in breast and gastric cancers, respectively), as measured by targeted MS, correlate with standard measures of HER2 positivity [1,2]. HER2 levels ≥ 2200 amol/ug were associated with longer OS [2].
- We hypothesized that targeted proteomics would identify high expression of HER2 (≥ 2200 amol/ug) in a range of tumor types, providing the potential to enroll in anti-HER2 clinical trials and may correlate with clinical benefit.

Methods





Figure 1: Microdissected tumor areas from FFPE tissue blocks (n=3,444) representing multiple indications were subjected to Liquid Tissue[®] and MS to quantify multiple protein targets in each patient sample.

Multiplexed Clinical Proteomics Menu

AR, ALK, AXL, CK-5, CK-7, EGFR, ERCC1, FGFR2, FR-alpha, hENT1, HER2, HER3, IGF1R, MET, MGMT, MSLN, PD-L1, p16, RON, ROS1, RRM1, SPARC, TOPO1, TOPO2A, TP63, TUBB3, TTF1

Results

HER2 expression above 2200 amol/ug in breast cancer correlates with increased survival



*Hazard ratio for OS cannot be determined because all patients with >2200 amol/ug of HER2 were alive after 6 years of anti-HER2 therapy

Figure 2. HER2 levels ≥ 2200 amol/ug were associated with longer OS (2).

HER2 super-expression in multiple tumor indications




Figure 3. HER2 ≥ 2200 amol/ug in multiple indications may be predictive of enhanced response to trastuzumab.

Response to trastuzumab in gallbladder cancer patient




Figure 4. A patient with gallbladder cancer had very high HER2 expression (3105 amol/ug) that responded to trastuzumab therapy.

Multiplex proteomic analysis identifies multiple therapeutic options in uterine cancer patient





Figure 5. A, B. Multiplex-MS analysis in uterine carcinoma patient identified multiple potential therapies.

Conclusions

- Personalized medicine should include the use of proteomic as well as genomic profiling to characterize all potential biomarkers which may direct therapies.
- Targeted Mass Spec-based proteomics offers a multiplexed, extremely sensitive, uniquely precise and highly reproducible method for tumor biomarker profiling and quantification to identify potential therapies.
- Patients with high HER2 protein expression (≥ 2200 amol/ug) in multiple cancer types benefit from anti-HER2 targeted therapies.
- Based on these and additional case studies, broader proteomic testing of tumor HER2 expression could extend the benefit of HER2 targeted therapy.

Reference

- Catenacci, D.V., et al.; Gastric Cancer, 2015.
- Nuciforo, P., et al.; Mol Oncol, 2016. **10**(1): p. 138-47.